



CALCIUM


O-cresolphthalein

Colorimetric test v/v

For in vitro diagnostic use only
Store at 2-8°C

PRINCIPLE OF THE METHOD

Calcium with o-cresolphthalein complex compound, at alkaline pH, yields a red colored complex, whose intensity is proportional to the calcium concentration.



CLINICAL SIGNIFICANCE

Calcium is the most abundant and one of the most important minerals in the human body. Approximately 99% of body calcium is found in bones. A decrease in albumin level causes a decrease in serum calcium. Among causes of hypercalcemia are cancers, large intake of vitamin D, enhanced renal retention, osteoporosis, sarcoidosis, thyrotoxicosis, hyperparathyroidism. Low levels of calcium are found in hypoparathyroidism, pseudohypoparathyroidism, vitamin D deficiency, malnutrition and intestinal malabsorption. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

| | | |
|---------------|--|--------------------------|
| R 1 Buffer | Ethanolamine | 500 mmol/L |
| R 2 Chromogen | o-Cresolphthalein 8-Hidroxyquinolein | 0.62 mmol/L 69 mmol/L |
| CALCIUM CAL | Calcium aqueous primary standard 10mg/dL | |

PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the reagents are ready for use. Stored at 2-8°C, it is stable up to the date of expiration as specified.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 570 nm ≥ 0.2 .

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 570 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

SAMPLES

- Serum or plasma: Separated from cells as rapidly as possible. Blood anticoagulants with oxalate or EDTA are not acceptable since these chemicals will strongly chelate calcium.
- Urine: Collect 24 hour urine specimen in calcium free containers. The collecting bottles should contain 10 ml of diluted Nitric acid (50% v/v). Record the volume. Dilute a sample 1/2 in distilled water. Mix. Multiply results by 2 (dilution factor). Stability of the samples: Calcium is stable 10 days at 2-8°C.

PROCEDURE

- Assay conditions:
Wavelength: 570 nm (550-590)
Cuvette : 1 cm. light path
Temperature 37°C / 15-25°C

- Adjust the instrument to zero with distilled water.
- Pipette into a cuvette:

| | Blank | Calibrator | Sample |
|------------------------------|-------|------------|--------|
| R 1 (mL) | 1.0 | 1.0 | 1.0 |
| R 2 (mL) | 1.0 | 1.0 | 1.0 |
| Calibrator (μL) | | 20 | -- |
| Sample (A) (μL) | -- | -- | 20 |

- Mix and incubate for 5 min. at 37°C / 15-25°C.
- Read the absorbance (A) of the samples and calibrator, against the Blank. The color is stable for at least 40 minutes.

CALCULATIONS

Serum and plasma

$$\frac{(A)\text{Sample}}{(A)\text{Calibrator}} \times 10 (\text{Calibrator conc.}) = \text{mg/dL}$$

calcium in the sample

Urine 24 h

$$\frac{(A)\text{Sample}}{(A)\text{Calibrator}} \times 10 (\text{Calibrator conc.}) \times \text{vol. (dL)}$$

urine/24 h = mg/24 h calcium

Conversion factor: mg/dL x 0.25 = mmol/L.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES' Serum or plasma:

Adults 8.5-10.5 mg /dL = 2.1-2.6 mmol/L
Children 10 -12 mg/dL = 2.5 - 3 mmol/L
Newborns 8 -13 mg/dL s 2 - 3.25 mmol/L
Urine:
Adults 50 - 300 mg/24h = 1.25 - 7.5 mmol/24h
Children 80 -160 mg/24h = 2 - 4 mmol/24h
These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 0.10 mg/dL to *linearity limit* of 15 mg/dL. If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

| | Inter-assay (n=20) | | Inter-assay (n=20) | |
|-------------|--------------------|------|--------------------|-------|
| Mean(mg/dL) | 9.08 | 15.7 | 9.03 | 14.28 |
| SD | 0.17 | 0.24 | 0.17 | 0.23 |
| CV% | 1.97 | 1.53 | 1.99 | 1.62 |

Sensitivity: 1 mg/dL = 0,034 A.

Accuracy Results obtained using ATLAS reagents (y) did not show systematic differences when compared with other commercial reagents (x). The results obtained using 50 samples were the following:

Correlation coefficient (r): 0.91.

Regression equation: $y=0.9069x + 0.9114$.

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

No interferences were observed with triglycerides up to 1.25 g/L.

A list of drugs and other interfering substances with calcium determination has been reported by Young et. al"

NOTES

- It is recommended to use disposable material. If glassware is used the material should be scrupulously cleaned with diluted 1/1 HNO₃ in water and then thoroughly rinsed it with distilled water.

2. Most of the detergents and water softening products used in the laboratories contains chelating agents. A defective rinsing will invalidate the procedure.
3. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
4. Use clean disposable pipette tips for its dispensation.

BIBLIOGRAPHY

1. Farell E C. Calcium. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1051-1255 and 418.
2. Kessler G. et al. Clin Chem Vol 10, No 8 1964; 686-706.
3. Connerty H. V. et al. Am J Clin Path Vol 45, No 3 1996; 200-296.
4. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
5. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACC 2001.
6. Burtis A. et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACC 1999.
7. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACC 1995.

ATLAS MEDICAL
Unit 4, William James House
Cowley Rd, Cambridge, CB4 0WX
Tel: ++44 (0) 1223 858 910
Fax: ++44 (0) 1223 858 524

PPI206A01
Revision A (05.07.2005)